



Surface disinfectant Metasys SD, refill 2

Surface disinfectant Metasys SD, refill 2, 1 x 5 litre canister

Available

SKU: 60030003

Categories: Surface Disinfection

Description

Surface disinfectant Metasys SD refill 2, 1 x 5 litre canister

Alcohol-based spray disinfection for alcohol-insensitive surfaces of medical devices.

- short action time;
- ready-to-use disinfection solution;
- formaldehyde-free;
- CE 0044.

Modern active agent combination: Ethanol, 2-propanol, quaternary ammonium compounds, glyoxal.