



## Surface disinfectant Metasys SD, refill 2

Surface disinfectant Metasys SD, refill 2, 1 x 5 litre canister

**Available**

**SKU:** 60030003

**Categories:** Surface Disinfection

### Description

---

Surface disinfectant Metasys SD refill 2, 1 x 5 litre canister

Alcohol-based spray disinfection for alcohol-insensitive surfaces of medical devices.

- short action time;
- ready-to-use disinfection solution;
- formaldehyde-free;
- CE 0044.

Modern active agent combination: Ethanol, 2-propanol, quaternary ammonium compounds, glyoxal.