



Surface disinfectant Metasys WD, start

Metasys WD surface disinfectant, start, 1 dispenser with 120 wipes.

Available on request

SKU: 60030018

Categories: Surface Disinfection

Description

Metasys WD surface disinfectant start, 1 dispenser with 120 wipes.

Alcohol-based wipes made of tear-resistant fabric, perfect for disinfecting small external surfaces and medical equipment.

Key Features:

- Sealed packaging: ensures product preservation and maximum hygiene.
- Long-lasting effect: prevents the settlement of new airborne germs.
- Advanced formula: a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

Broad-spectrum action:

- HBV, HIV, HCV: disinfects in 30 seconds.
- Bactericidal and fungicidal: effective in 1 minute.
- Tuberculocidal: action in 5 minutes.

Certification: CE 0044.

Packaging Metasys WD surface disinfectant start: 1 dispenser with 120 wipes.