



Surface disinfectant Metasys WD, start

Metasys WD surface disinfectant, start, 1 dispenser with 120 wipes.

Shipping expected in approximately 15 days

SKU: 60030018

Categories: Surface Disinfection

Description

Metasys WD surface disinfectant start, 1 dispenser with 120 wipes.

Alcohol-based wipes made of tear-resistant fabric, perfect for disinfecting small external surfaces and medical equipment.

Key Features:

- **Sealed packaging:** ensures product preservation and maximum hygiene.
- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** a combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.

Broad-spectrum action:

- **HBV, HIV, HCV:** disinfects in 30 seconds.
- **Bactericidal and fungicidal:** effective in 1 minute.
- **Tuberculocidal:** action in 5 minutes.

Certification: CE 0044.

Packaging Metasys WD surface disinfectant start: 1 dispenser with 120 wipes.