



## Dental cart Prophy 900, 3 instruments

Dental cart Prophy 900, 3 instruments: 1x Luzzani Minimate steel syringe, 1x LED induction micromotor, 4-way turbine cord, safety handpiece.

**Available on request**

**SKU:** DE900E

**Categories:**

### Description

---

Prophy 900 surgical cart has all the tools and features you need for surgery and is designed to ensure maximum operational efficiency.

The setup of the available tools can be customized based on the specific clinical needs: Prophy 900 features a instrument table that can be set up for maximum five different instruments.

The quality of the materials used ensures utmost sturdiness and cutting-edge design: the 2cm thick polyurethane instrument table is covered with an epoxy paint resistant to disinfectants. Unlike plastics, the Prophy 900 coating does not turn yellow and does not crack.

- Maximum efficiency: it guarantees smooth work processes, providing all the necessary tools.
- Designed to ensure operational efficiency: it is equipped with all the functions necessary to carry out the interventions.
- All parameters can be perfectly kept under control thanks to the wide LCD control panel with optional soft touch keys and progressive control pedal.
- Indestructible: 2 cm thick polyurethane instrument table.
- Easy movement thanks to the comfortable anti-static casters and the pedal shelf

The separate water supply guarantees total operating soprattutto in ambito chirurgico. autonomy; there is also an independent spray adjustment for each instrument.

Thanks to the 5 compartments of the instrument table, Prophy 900 surgical cart can be set up with the most suitable equipment for every need.

Prophy 900 features a large table top on which are inserted two tool trays, 175x280x15 mm, removable and autoclavable in order to ensure maximum hygiene, which is absolutely necessary for surgery.