



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/surgical-green-contra-angle-ivory-silk-line-201-external-spray/>



Surgical green contra-angle Ivory Silk Line 20:1 external spray

Surgical green contra-angle Ivory Silk Line 20:1, external spray

400,80 €

Available

SKU: SL4201C

Categories: Professional Equipment

Description

The Surgical green contra-angle Ivory Silk Line 20:1, external spray, has been developed to respond to the specific demands of oral, stomatological and facial surgery, as well as for implantology.

Without Fiber optic, thanks to the special structure of noble metal placed in the front that increases the strength of the instrument to preserve it from the continuous tensions and to its rear part that balances the weight lightening it perfectly, the SL4201C with external spray represents the solution that combines solidity and ergonomics.

In addition, the design of the SL4201C consists of two parts that make it easier to disassemble and thus facilitate cleaning and maintenance, while maintaining its aesthetic and ergonomic properties.

Technical features Surgical green contra-angle Ivory Silk Line 20:1

Removable head	NO
Optical fiber	NO
Transmission	20:1
Rotor bearings	In Ceramics
Indications	Implantology
Stop mechanism	Button
Type of Spray	External
Compatibility	Intramatic Universal
Rotation speed	2.000
Head Rotation	NO
Torque (Ncm):	? 55
Diameter of the Head	9.5mm
Head Height	13.5mm
Diameter of the Bur	2.334/2.350 mm

Length of the Bur	21.95mm
Weight	58gr
Autoclavable	YES
Type of Lubrication	Oil