



# Surgical green contra-angle Ivory Silk Line 20:1 external spray

Surgical green contra-angle Ivory Silk Line 20:1, external spray

**Available** 

**SKU: SL4201C** 

Categories: Contra-Angles

# Description

# Ivory Silk Line Green Surgical Contra-Angle 20:1 Blind, External Spray

The **Ivory Silk Line green surgical contra-angle 20:1 blind external spray** is designed to meet the specific needs of oral, maxillofacial and implant surgery. Its precision microtechnology ensures excellent handling, reliability and full control in every procedure.

## High-efficiency mechanical transmission.

The internal mechanical system is built with **extremely low tolerances**, transferring micromotor power efficiently without energy loss or excess noise. This eliminates vibration and provides smooth, precise and stable performance.

## Compact head for optimal access.

The **small head design** allows easy access to narrow and sensitive areas, offering excellent control and precision even in complex surgical conditions.

### Sterilisable and durable design.

All Ivory Silk Line contra-angles are **thermodisinfectable and autoclavable up to 135°C**, ensuring total hygiene, long-term reliability and consistent clinical performance.

#### Technical features

Removable head	NO
Optical fiber	NO



Transmission	20:1
Rotor bearings	In Ceramics
Indications	Implantology
Stop mechanism	Button
Type of Spray	External
Compatibility	Intramatic Universal
Rotation speed	2.000
Head Rotation	NO
Torque (Ncm):	? 55
Diameter of the Head	9.5mm
Head Height	13.5mm
Diameter of the Bur	2.334/2.350 mm
Length of the Bur	21.95mm
Weight	58gr
Autoclavable	YES
Type of Lubrication	Oil