



## Surgical green contra-angle Ivory Silk Line 20:1 external spray

Surgical green contra-angle Ivory Silk Line 20:1, external spray

**Available**

**SKU:** SL4201C

**Categories:** Contra-Angles handpieces

### Description

#### Ivory Silk Line Green Surgical Contra-Angle 20:1 Blind, External Spray

The **Ivory Silk Line green surgical contra-angle 20:1 blind external spray** is designed to meet the specific needs of oral, maxillofacial and implant surgery. Its precision microtechnology ensures excellent handling, reliability and full control in every procedure.

#### High-efficiency mechanical transmission.

The internal mechanical system is built with **extremely low tolerances**, transferring micromotor power efficiently without energy loss or excess noise. This eliminates vibration and provides smooth, precise and stable performance.

#### Compact head for optimal access.

The **small head design** allows easy access to narrow and sensitive areas, offering excellent control and precision even in complex surgical conditions.

#### Sterilisable and durable design.

All Ivory Silk Line contra-angles are **thermodisinfected and autoclavable up to 135°C**, ensuring total hygiene, long-term reliability and consistent clinical performance.

#### Technical features

Removable head	No
----------------	----



Fiber optics	No
Transmission ratio	20:1
Indications	Implantology
Locking mechanism	Push-button
Spray type	External
Compatibility	Intramatic ISO 3964
Rotation speed	2000rpm
Torque	? 55Ncm
Rotatable head	No
Head diameter	9.5mm
Head height	13.5mm
Bur diameter	2.334/2.350mm
Bur length	9–34mm
Weight	58g
Autoclavable	Yes (max 135°C)
Lubrication	Oil