



Surgical green contra-angle Ivory Silk Line 20:1 external spray

Surgical green contra-angle Ivory Silk Line 20:1, external spray

Shipping expected in approximately 15 days

SKU: SL4201C

Categories: Contra-Angles handpieces

Description

The **Ivory Silk Line 20:1 surgical contra-angle with external spray** is designed to support implantology procedures in the dental practice, offering control, precision and stability during clinical work.

The **20:1 reduction transmission** allows operation at a controlled speed with high torque, making the instrument suitable for surgical procedures that require precise movement control and reliable performance.

The high-precision mechanical transmission components are manufactured with extremely reduced tolerances, supporting efficient power transfer from the micromotor. This helps reduce vibrations, noise and energy loss, providing smoother, more stable and more precise operation during procedures.

The **compact head** facilitates access to surgical areas that are more difficult to reach, improving visibility and control during treatment.

External spray system contributes to cooling the working area during use, supporting operating safety in implantology procedures.

The contra-angle is **thermo-disinfectable and sterilizable in an autoclave up to 135°C**, in accordance with the hygiene protocols required in dental practices. This feature ensures complete hygiene, long service life and consistent safety, even during intensive use.

Technical features Ivory Silk Line 20:1 surgical contra-angle

Removable head	No
Fiber optics	No
Transmission ratio	20:1
Indications	Implantology



Locking mechanism	Push-button
Spray type	External
Compatibility	Intramatic ISO 3964
Rotation speed	2000rpm
Torque	? 55Ncm
Rotatable head	No
Head diameter	9.5mm
Head height	13.5mm
Bur diameter	2.334/2.350mm
Bur length	9–34mm
Weight	58g
Autoclavable	Yes (max 135°C)
Lubrication	Oil