



## Surgical multiLED 2 lamps, 160.000 +160.000 Lux, ceiling

Surgical multiLED 2 operating lamp, dual LED head 160,000 + 160,000 Lux, ceiling mounted

**Shipping expected in approximately 60 days**

**SKU:** DE6.163

**Categories:** Surgical Lights

### Description

The Surgical multiLED 2 operating lamp 160.000 + 160.000 Lux ceiling (DOUBLE HEAD) offer excellent illumination for all dental and non-dental surgical needs.

Functional and linear design, combined with a lightweight structure, make for unrivalled easy-to-handle lights.

The suspension system with balanced articulated arm with 6 joints allows stable and precise positioning.

All SURGICAL multiLED models are equipped with a colour temperature compensation technology ranging from 3300 K° to 5000 K° and manual adjustment of the light intensity using the control panel.

They are also equipped with the “Endo mode” function; that is the support light for endoscopy operations offering light intensity and colour temperature selection.

Light intensity can easily be seen on the v-meter bar located on the control panel. The total reflection optical system is made with CAD/CAM technology, its illumination depth is up to 1300 mm.

### Technical features Surgical multiLED 2 operating lamp 160.000 + 160.000 Lux ceiling (DOUBLE HEAD)

Illuminance (lux)	160.000 + 160.000
Colour temperature (k°)	from 3300 to 5000
Colour rendering index (Ra)	>95
Illumination depth (mm)	?1300



Total irradiance (W/m2)	426
Bright field size (mm)	160 ~ 280
LED bulb lifespan (h)	50,000
Power supply (V/Hz)	AC110-240v, 50/60Hz
Total energy absorption (W )	140 W
Maximum installation height (metres)	2,9
Certification	2007/47/CE – IEC62321
Compliance with	(UE) 2017/745 standards