



## Antibacterial medical keyboard, autoclavable. IP65

Antibacterial medical keyboard, IP65, with autoclavable silicone key membrane, protects the key field against ingress of dirt and fluids.

**Estimated shipping in 10 days**

**SKU:** DE1004.5SO

**Categories:** All-in-One Medical PCs

### Description

---

The Autoclavable antibacterial medical keyboard IP65 is a high quality, highly reliable notebook PC keyboard.

The purification and disinfection capacity makes it a perfect fit in hospitals and at practitioners, microbiological laboratories as well as industrial environments.

The silicone key membrane protects the key field against ingress of dirt and fluids.

The Antibacterial medical keyboard impresses with precise key actuation and accurate key guidance.

It is highly recommendable for heavy duty type writing in dirt and hygiene critical operations.

#### **DECISIVE BENEFITS Autoclavable antibacterial medical keyboard IP65**

- Thorough wipe disinfection of key field and keyboard
- Silicone key membrane with low key profile and self-containing, seamless surface
- CleanFunction allows to disable key field while cleaning ("Fn"+"Clean on/off")
- Exchangeable key membrane for separate reprocessing (Disinfection, Steam Sterilization)
- Suitable for a multitude of hospital cleaning agents and surface disinfectants (VAH)
- Protection IP65