



Tecnomed Digiclean stainless steel ultrasound tank, 10 l with lid

Tecnomed Digiclean stainless steel ultrasonic tank, 10 l capacity, digital timer, with lid and drain cock

Shipping expected in approximately 30 days

SKU: CLEA10L

Categories: Ultrasonic Cleaners

Description

Ultrasonic dental cleaner Digiclean 10 litres.

The Ultrasonic washing is the most technologically advanced and cost-effective means of cleaning.

Ultrasonic cleaner Digiclean 10l generate an high ultrasonic power that associated with a disinfectant cleanser will ensure you a quick and effective cleaning action by removing dirt even in the hard- to-access parts.

Made of non-scratch stainless steel, it proves resistant over time against wear-and-tear corrosive agents.

An electronic generator operated by a microprocessor optimizes power performance while a display control panel allows the setting and the viewing of vibration times and liquid temperature values.

- An Intuitive control panel lets you set cleaning time (from 1 to 99 minuts), solution temperature (0°C to 80°C) and ultrasonic frequency.
- Our ultrasonic cleaners are made entirely from stainless steel for maximum durability.
- The 5 liters and 10 liters models are equipped with an exhaust washing liquid faucet so you can drain water quickly and easily.
- A new feature of this range of cleaners is the availability of Pyrex beakers (optional) for indirect cleaning with particularly aggressive solutions.

Comparative table Ultrasonic dental cleaner Digiclean 10 litres.

FEATURES	Code CLEA3L	Code CLEA5L	Code CLEA10L



nominal capacity	3,04 liters (237x93x138)mm	5,04 liters (240x150x140)mm	10,8 liters (300x150x240)mm
Real Capacity	3 liters	4,350 liters	9,5 liters
Transducers	2	2	4
Swinging power	120	120	240
Swinging frequency	40KHZ	40KHZ	40KHZ
Drain tap	?	?	?
Digital timer	?	?	?
Heating	?	?	?
Large/medium basket	?	?	?
Small basket	?	?	?
Cover	?	?	?
Cups support	Opt.	Opt.	Opt.
Cups	Opt.	Opt.	Opt.
Power	230 vac	230 vac	230 vac
Weight gr.	4 Kg	5 Kg	10 Kg



External measures lxdxl mm	252 x 150 x 271	252 x 150 x 271	312 x 250 x 337
Warranty	24 months	24 months	24 months