



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/tecnomed-digiclean-stainless-steel-ultrasonic-tank-5l-with-lid/>



Tecnomed Digiclean stainless steel ultrasonic tank, 5 l with lid

Tecnomed Digiclean stainless steel ultrasonic tank, 5 l capacity, digital timer, with lid and drain cock

591,16 €

Available

SKU: CLEA5L

Categories: Ultrasonic cleaners

Description

Ultrasonic cleaner Digiclean 5l.

The Ultrasonic washing is the most technologically advanced and cost-effective means of cleaning.

Ultrasonic cleaner Digiclean 5l generate an high ultrasonic power that associated with a disinfectant cleanser will ensure you a quick and effective cleaning action by removing dirt even in the hard- to-access parts.

Made of non-scratch stainless steel, it proves resistant over time against wear-and-tear corrosive agents.

An electronic generator operated by a microprocessor optimizes power performance while a display control panel allows the setting and the viewing of vibration times and liquid temperature values.

- An Intuitive control panel lets you set cleaning time (from 1 to 99 minuts), solution temperature (0°C to 80°C) and ultrasonic frequency.
- Our ultrasonic cleaners are made entirely from stainless steel for maximum durability.
- The 5 liters and 10 liters models are equipped with an exhaust washing liquid faucet so you can drain water quickly and easily.
- A new feature of this range of cleaners is the availability of Pyrex beakers (optional) for indirect cleaning with particularly aggressive solutions.

Comparative table Ultrasonic cleaners Digiclean:

FEATURES	Code CLEA3L	Code CLEA5L	Code CLEA10L
nominal capacity	3,04 liters (237x93x138)mm	5,04 liters (240x150x140)mm	10,8 liters (300x150x240)mm

Real Capacity	3 liters	4,350 liters	9,5 liters
Transducers	2	2	4
Swinging power	120	120	240
Swinging frequency	40KHZ	40KHZ	40KHZ
Drain tap	?	?	?
Digital timer	?	?	?
Heating	?	?	?
Large/medium basket	?	?	?
Small basket	?	?	?
Cover	?	?	?
Cups support	Opt.	Opt.	Opt.
Cups	Opt.	Opt.	Opt.
Power	230 vac	230 vac	230 vac
Weight gr.	4 Kg	5 Kg	10 Kg
External measures lxdxh mm	252 x 150 x 271	252 x 150 x 271	312 x 250 x 337
Warranty	24 months	24 months	24 months