



Tecnomed Digiclean ultrasound tank 5 lt, REFURBISHED 2 YEARS WARRANTY

Tecnomed Digiclean stainless steel ultrasonic tank, 5 l capacity, digital timer, with lid and drain cock, refurbished product with 2 years warranty.

Available

SKU: OFOUT128

Categories: Outlet

Description

Ultrasonic cleaner Digiclean 5l. REFURBISHED PRODUCT, 2 YEARS WARRANTY.

The Ultrasonic washing is the most technologically advanced and cost-effective means of cleaning.

Ultrasonic cleaner Digiclean 5l generate an high ultrasonic power that associated with a disinfectant cleanser will ensure you a quick and effective cleaning action by removing dirt even in the hard- to-access parts.

Made of non-scratch stainless steel, it proves resistant over time against wear-and-tear corrosive agents.

An electronic generator operated by a microprocessor optimizes power performance while a display control panel allows the setting and the viewing of vibration times and liquid temperature values.

- An Intuitive control panel lets you set cleaning time (from 1 to 99 minuts), solution temperature (0°C to 80°C) and ultrasonic frequency.
- Our ultrasonic cleaners are made entirely from stainless steel for maximum durability.
- The 5 liters and 10 liters models are equipped with an exhaust washing liquid faucet so you can drain water quickly and easily.
- A new feature of this range of cleaners is the availability of Pyrex beakers (optional) for indirect cleaning with particularly aggressive solutions.

Comparative table Ultrasonic cleaners Digiclean:

Features	CLEA5L
Nominal capacity	5,04 liters (240 x 150 x 140) mm
Real capacity	4,350 liters
Transducers	2



Features	CLEA5L
Swinging power	120 W
Swinging frequency	40 KHz
Drain tap	?
Digital timer	?
Heating	?
Large/medium basket	?
Small basket	?
Cover	?
Cups support	Opt.
Cups	Opt.
Power	230 Vac
Weight	5 Kg
External measures (l x d x h)	252 x 150 x 271 mm
Warranty	24 months