



Termometro a infrarossi reg.Ministero Salute. DISPONIBILITÀ: 7 pz

Professional infrared thermometer for fast, accurate and hygienic body temperature measurement without direct contact.

Available

SKU: DE2003

Categories: Professional Equipment

Description

The **TERMOTEC PRO infrared thermometer** is a professional non-contact device designed for the fast and accurate measurement of body temperature before examinations, treatments and clinical procedures. Developed for medical environments, it allows temperature readings to be taken hygienically and efficiently without direct contact with the patient.

Suitable for use at a distance of up to 5cm, it provides reliable temperature measurement within seconds. The backlit display changes colour according to the recorded value, making readings clear and immediate while helping operators identify normal, elevated or high temperatures at a glance.

The device also features an alarm system to indicate any errors during temperature acquisition, along with an audible alert when body temperature exceeds 37.5°C. Its ergonomic shape and intuitive controls make it practical and easy to use in busy professional settings.

Technical features

- Professional non-contact infrared thermometer
- Maximum measuring distance: 5cm
- Automatic shut-off
- Dimensions: 145x88x42mm
- Weight: 140g
- Body temperature operating range: 32.1 to 42.9°C
- Power supply: 2 AAA batteries, DC 3.0V
- Certification: 93/42EEC
- Class: IEC60601