


Titanium scaler tip, perio P56T Woodpecker compatible



Perio tip for scaler, for periodontal treatments, Woodpecker P56T compatible.

Available

SKU: PG056T

Categories: EMS - Woodpecker

Description

The titanium perio tip compatible with Woodpecker P56T is a scaler tip designed for professional periodontal treatments.

Its slim shape allows access to deep periodontal pockets and supports treatment in difficult-to-reach areas, helping remove biofilm and deposits along the root surface.

Specifications

- Compatible with Woodpecker P56T
- Suitable for use with Woodpecker-compatible piezoelectric scalers
- For professional dental use