


Titanium scaler tip, perio P56T Woodpecker compatible



Perio tip for scaler, for periodontal treatments, Woodpecker P56T compatible.

Available

SKU: PG056T

Categories: EMS - Woodpecker

Description

Titanium scaler tip perio, compatible Woodpecker P56T.

Tecnomed Italia inserts for periodontal treatments are extremely thin and able to enter even deep pockets easily, enabling you to achieve excellent and painless results in radicular debridement.

They are also ideal for removing tartar from tight interdental spaces and touching crowns. The reduced oscillations generated by the low power settings needed for these tips makes treatment with them particularly pleasant for the patient.

Their subgingival spray system eliminates biofilm fragments more effectively.

Compatibility Titanium scaler tip perio: Woodpecker P56T.