



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/titanium-scaler-tip-perio-p59t-woodpecker-compatible/>

## Titanium scaler tip, perio P59T Woodpecker compatible


Perio tip for scaler, for periodontal treatments,  
Woodpecker P59T compatible.

**55,58 €**

**Available**

**SKU:** PG059T

**Categories:** Periodontal inserts



### Description

Titanium scaler tip perio PG059T.

Tecnomed Italia inserts for periodontal treatments are extremely thin and able to enter even deep pockets easily, enabling you to achieve excellent and painless results in radicular debridement.

They are also ideal for removing tartar from tight interdental spaces and touching crowns.

The reduced oscillations generated by the low power settings needed for these tips makes treatment with them particularly pleasant for the patient.

Their subgingival spray system eliminates biofilm fragments more effectively.

Compatibility Titanium scaler tip perio PG059T.

- Woodpecker P59T.