



TMI laboratory micromotor, universal brush drive, 50.000 rpm



Tecnomed Italia laboratory micromotor. Max. speed 50,000rpm, induction. Smooth cord.

Available

SKU: MX950

Categories: Handpieces

Description

TMI laboratory micromotor 50.000 rpm.

Tecnomed's micromotor handpieces are ergonomically designed and offer a perfect balance, not only of the handpiece alone but of the complete handpiece and cord system.

Because they do not pull on your wrist, your working position can be more relaxed and natural and you can concentrate better on the intricate, complex and time-consuming tasks you often have to perform.

Maintenance is effortless too, thanks to an easy to clean spindle and quick to replace bearings. Excellent speed control, mechanical robustness and ergonomic design are just some of the key features of these handpieces.

They are also fully interchangeable and can be used with any dental unit that provides the right output voltage. The choice is yours!

Compatible with brushless motor control units*

*In some cases it may be necessary to change the micromotor handpiece connector.

Technical features TMI laboratory micromotor 50.000 rpm.

| | |
|----------------------|------------|
| Micromotor handpiece | brushless |
| Maximum speed | 50.000 rpm |
| Handpiece lenght | 162 mm |
| Handpiece diameter | 26 mm |



| | |
|------------------|---------|
| handpiece weight | 358 gr. |
| Bur diameter | 2,35 mm |