



Printed on 30/05/2025

<https://shop.tecnomeditalia.com/en/product/universal-pyrex-glass-beaker-for-ultrasonic-baths-600-ml/>



Universal Pyrex glass beaker for ultrasonic baths, 600 ml

Universal Pyrex glass beaker for ultrasonic baths, capacity 600 ml, with o-ring support

18,27 €

Available

SKU: CLA03

Categories: Ultrasonic cleaners

Description

Universal Pyrex glass beaker.

Universal Pyrex glass beaker for ultrasonic baths, 600 ml, equipped with o-ring for support.

Diameter 100mm, Height 125mm.