



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/wek-light-metasy-s-water-decontamination-system-free-standing/>



Wek Light Metasys water decontamination system, free-standing

Wek Light Metasys water decontamination system, without DVGW, floor-mounted

1.651,10 €

Available on request

SKU: 05020026

Categories: Water decontaminators

Description

Free-standing Wek Light Metasys water decontamination system is the ideal solution to maintain optimal hygiene in the water lines of dental units, safeguarding both patients and instruments from potential contamination.

Designed to prevent biofilm formation and ensure maximum cleanliness of the pipelines, the WEK Light eliminates the risk of bacterial growth and unpleasant odors, enhancing the safety and efficiency of dental treatments.

Key Features / Benefits:

- **Automatic dosage control:** Ensures precise delivery of the decontaminating solution, optimizing effectiveness without waste.
- **Maintenance-free:** Requires no complex maintenance; the only task needed is replacing the decontaminant solution cartridges when depleted.
- **Constant protection:** Keeps water lines in a state of maximum hygiene, preventing biofilm formation and extending the lifespan of your dental equipment.
- **Regulatory compliance:** Designed to meet current hygiene and safety standards, providing quality assurance for any dental practice.

Technical Specifications:

- **Certification:** CE 0044
- **Installation type:** Floor
- **Dosage system:** Automatic
- **Maintenance:** No maintenance required, only cartridge replacement
- **Compatibility:** Suitable for any dental unit
- **Dimensions:** Compact for easy integration into the practice.

The WEK Light Water Decontaminator by Metasys floor model is your reliable partner for a safer and compliant work environment.

Choose the WEK Light products by Metasys from Tecnomed Italia and discover a new level of efficiency and reliability in water decontamination systems for dental practices. Optimize your daily work with the most advanced suction technology on the market, with the assurance and guarantee of Tecnomed Italia's after-sales service.

Technical features Free-standing Light Wek Metasys water decontamination system.

Model	Wek light
Installation	Own cabinet
Product code for hydrogen peroxide solution	60040100-1
Power supply	24VAC
Power	180 Watt
Weight	6.3 kg
Dimensions (w x l x h) mm	200 x 155 x 315 mm
Warranty	12 months