



STERILA surgical sink, white stainless steel with UV-C disinfection, 800

White surgical sink with UV-C sterilization and hands-free electronic control for clinical environments. Dimensions: 800 x 650 x 1020 mm

Available on request

SKU: WASHCH-380B

Categories: Sinks

Description

White surgical sink with UV-C sterilization and electronic control

The white surgical sink is designed to ensure high hygiene standards in clinical and surgical environments, while offering improved aesthetic integration within dental practices and modern healthcare facilities.

The stainless steel structure with white coated finish allows seamless integration with medical furniture and dental units, while maintaining the technical performance required for high hygiene environments.

UV-C sterilization for advanced microbiological control

The integrated UV-C sterilization system is designed to reduce microbial load within the water circuit, improving hygiene safety directly at the point of use. This makes it particularly suitable for surgical and clinical environments.

Hands-free electronic control

The sink is equipped with an electronic system with solenoid valves, allowing water activation without direct contact. This hands-free operation reduces cross-contamination risk and improves operational efficiency.

Retractable design for space optimisation

The retractable configuration allows the sink to be closed or concealed when not in use, reducing visual clutter and optimising workspace organisation. This is especially beneficial in dental practices and modern clinical environments.



White finish for clinical integration

The white coated finish enhances visual cleanliness and allows perfect integration with dental units and clinical furniture, contributing to a professional and organised environment.

This feature is particularly relevant in environments where aesthetics influence patient perception and overall experience.

Medical-grade materials

- corrosion-resistant structure
- compatible with medical disinfectants
- easy cleaning and sanitisation
- long-lasting performance

Applications

- dental practices
- medical clinics
- outpatient facilities
- operating rooms
- high hygiene environments

Operational advantages

- integrated UV-C sterilization
- hands-free electronic control
- reduced contamination risk
- space-saving retractable design
- improved aesthetic integration
- easy-to-clean surfaces

Technical specifications

Type: retractable surgical sink

Material: stainless steel with white coated finish

Control system: electronic with solenoid valves

Sterilization: integrated UV-C system

Installation: clinical and surgical environments

Dimensions: 800 × 650 × 1020 mm (H)