



VetLine Prophy 600 4-hole veterinary air polishing handpiece

VetLine Prophy 600 veterinary air polishing handpiece with 4-hole connection, designed for professional veterinary prophylaxis, plaque removal and animal dental cleaning.

Available

SKU: PL600H

Categories: Veterinary Handpieces

Description

The **VetLine Prophy 600 4-hole veterinary air polishing handpiece** is a professional sodium bicarbonate powder device designed to support **veterinary prophylaxis and dental hygiene procedures** in practices and clinics.

Intended for **veterinary use**, Prophy 600 helps remove plaque, superficial staining and deposits from the tooth surface, making it suitable for animal dental cleaning as part of routine oral hygiene procedures.

The nozzle design helps direct the powder, air and water jet in a controlled way, reducing powder dispersion and supporting continuous rinsing of the oral cavity during treatment.

The **360° rotating tip** allows easier access to difficult-to-reach areas, helping the operator perform accurate cleaning even in less accessible zones of the mouth.

The tip is suitable for **thermal disinfection and autoclaving at 135°C**, supporting high hygiene standards during professional veterinary use.

TECHNICAL FEATURES

Product type: veterinary air polishing handpiece

Line: VetLine

Model: Prophy 600

Intended use: veterinary use

Connection: 4-hole

Air supply: 2.8 bar

Water supply: 1-2 bar

Working system: sodium bicarbonate powder, air and water

Tip: 360° rotating tip

Tip sterilisation: autoclavable at 135°C



Use: veterinary prophylaxis, plaque removal and animal dental cleaning

Replacement O-ring kit: included