



Wireless electric ulv nebulizer foruse with GREEN&CLEAN SK, 800 ml

The Flower Clean 800ml wireless nebulizer is designed for fast environmental sanitisation, with an **18 Volt 4 Ah battery**, **800ml tank**, **adjustable 50-120 micron droplet size** and **output up to 350ml/min**.

Available

SKU: PC5007-3

Categories: Professional Equipment

Description

Flower Clean 800ml wireless nebulizer is a practical solution for the fast and even sanitisation of surfaces and indoor environments in professional medical settings. Designed for straightforward and immediate use, it delivers the sanitising product effectively, helping to reduce both treatment time and operating costs in routine disinfection procedures.

Thanks to its **800ml tank** and **adjustable droplet size from 50 to 120 microns**, the unit can be adapted to different application requirements, ensuring controlled and consistent liquid distribution. The **adjustable output up to 350ml/min** and **spraying distance up to 1.5m** make it suitable for treating rooms and surfaces quickly and efficiently.

This model is powered by an **18 Volt 4 Ah battery**, delivers **550W** of power and operates on a discontinuous cycle of **10 minutes on and 10 minutes off**. Its wireless design allows greater freedom of movement during use, making it especially practical in professional environments where mobility and ease of handling are important during sanitisation procedures.

It is an effective choice for operators looking for a reliable, easy-to-use nebulizer suitable for fast environmental sanitisation in medical and professional settings.

Benefits

- Reduces disinfection time and operating costs
- Easy and immediate to use
- Adjustable nebulisation according to application needs
- Good working autonomy thanks to the 800ml tank
- Greater freedom of movement thanks to the wireless design
- Effective coverage up to 1.5m

Technical specifications



Specification	Value
Product	Flower Clean 800ml wireless nebulizer
Power supply	18 Volt 4 Ah battery
Current	2.5A
Power	550W
Tank capacity	800ml
Droplet size	Adjustable 50-120 microns
Spraying distance	Up to 1.5m
Adjustable output	Up to 350ml/min
Dimensions with tank	290x275x115
Weight with full tank approx.	2.5kg
Power cable length	/
Operating mode	Discontinuous, 10 minutes on, 10 minutes off
Certification	CE
Warranty	24 months