

## Woodpecker PT-3 led ultrasonic scaler, EMS compatible



Woodpecker led benchtop power scaling unit with tank, 2 handpieces and 16 tips, circular vibrations, EMS compatible

**Shipping expected in approximately 30 days**

**SKU:** SCR3000

**Categories:** Periodontal Scalers

### Description

---

#### Woodpecker PT-3 led ultrasonic scaler Compatible with EMS

The PT-3 LED periodontal scaler by Woodpecker is the ideal choice for those seeking an innovative device designed to provide efficient periodontal treatments while minimizing patient discomfort and ensuring consistent high-quality results.

Equipped with ceramic piezoelectric ultrasound technology, the PT-3 scaler delivers high performance with optimal precision, reducing the risk of tooth damage compared to traditional instruments. The ergonomic handpiece with integrated LED ensures excellent visibility, allowing easy access to even the most challenging areas, making the practitioner's work more effective.

#### Key Features and Strengths:

- **Elliptical Movement:** With elliptical vibrations ranging between 20 and 60  $\mu\text{m}$ , the scaler offers a gentle, painless treatment for patients, reducing operational time and improving treatment effectiveness.
- **Titanium Inserts:** Specifically designed for delicate periodontal treatments, the titanium inserts ensure optimal vibration transmission while protecting enamel, offering maximum comfort for patients.
- **Smart Touch Display:** The intuitive display allows for precise control of all functions, and its water-resistant design makes cleaning and disinfection quick and easy after treatments.
- **“Easy & Fast” Refill:** The graduated tank is easy to refill and provides immediate visibility of the liquid level, optimizing the management of the device during procedures.

**Ideal for:** The PT-3 LED is particularly suited for dentists and periodontists who need a reliable, efficient tool for precise periodontal treatments. It ensures professional results with minimal discomfort, enhanced patient satisfaction, and improved treatment efficiency.



Choose PT-3 LED scaler by Woodpecker from Tecnomed Italia and discover a new level of efficiency and reliability in dental treatments. Optimize your daily work with DTE-Woodpecker solutions, with the assurance and guarantee of Tecnomed Italia's after-sales service.

#### **Technical Specifications Woodpecker PT-3 led ultrasonic scaler:**

- Technology: Ceramic piezoelectric ultrasound
- Compatibility: EMS
- Vibration amplitude: 20-60  $\mu\text{m}$
- Inserts: Titanium
- Handpiece: LED integrated
- Display: Smart touch
- Tank: Graduated with quick refill system